
CRYO PRO


- EN • INSTRUCTIONS FOR USE
- FR • MANUEL D'UTILISATION
- DE • BEDIENUNGSANLEITUNG
- IT • ISTRUZIONI PER L'USO
- ES • MANUAL DE EMPLEO
- NL • GEBRUIKERSHANDLEIDING

TABLE OF CONTENTS

(1) Composition	(9) Warranty
(2) Warning and cautions	(10) Disposal and recycling
(3) Intended use	(11) Troubleshooting
(4) Contraindications	(12) Specifications
(5) Operating the system	(13) Symbols used
(6) Cleaning and storage	(14) Information about electromagnetic compatibility (EMC)
(7) Cuffs	
(8) Accessories	

(1) COMPOSITION (DA03020001)

① Control Unit	
② Connector Hose	
③ AC Adaptor	
④ Cuff (not included) see § (7) Cuffs	

(2) WARNING AND CAUTIONS

READ THE INSTRUCTIONS BEFORE USE
TO AVOID ANY DANGER DUE TO MISUSE

- Do not apply the cuff directly to skin, add an isolation layer between the cuff and the skin.
- Only use Cryo cuffs, power adaptor, spare parts and accessories approved and supplied by the manufacturer.
- Do not drop or stack heavy items on the product.
- Do not modify, disassemble or destroy the product and only use for its intended purpose.
- Do not use the product when water leaks.
- Stop use immediately if the product or electrical parts are damaged / leaking.
- Be careful not to trip over the system's power cords or hose.
- To avoid damage, do not operate the system if water is not filled in the cooler.

- The control unit is not refrigeration equipment, need to add ice water or ice cubes to adjust the temperature.
- Keep the cooler away from sunlight and heat source.
- To avoid risk of electrical shock, always turn off the system and disconnect power cord from its electrical outlet when not in use or before adding or emptying ice and water.
- The product should only be used following the treatment recommendations of your healthcare practitioner.
- Children, the elderly and those who have cognitive disabilities or communication barriers cannot use the product alone without supervision.
- This product is forbidden to be used on any wound.
- When using the product, the host machine and the ice capsule should be placed at the same level, and the height difference between the two should not be more than 40cm.
- Monitor the skin condition of affected area for any burning, itching, increased swelling, or pain. If any of these happens, or any changes in skin appearance occur such as blisters, increased redness or discoloration, users should stop using and consult a physician.

(3) INTENDED USE

Dr Aktive® Cryo Pro pulse cold therapy system is intended for physical cooling, cold treatment and joint protection. It is intended for postoperative care and acute injury to reduce oedema, swelling and pain.

Dr Aktive® Cryo Pro is intended to be used by or on the order of licensed healthcare professionals used in hospitals, outpatient clinics, athletic training or home settings.


(4) CONTRAINDICATIONS

DO NOT USE THE PRODUCT BEFORE CONSULTING A DOCTOR
IF ANY OF THE FOLLOWING APPLIES


- People with Raynaud's disease.
- People with acute, unstable (untreated) fracture.
- People with rheumatoid arthritis in the affected region.
- People with open wounds, or localized unstable skin condition.
- People who are intolerant to cold compress or have cold hypersensitivity.
- The pregnant, people with hypertension, extreme low blood pressure or circulation disorders.
- People with vascular diseases such as prior frostbite, diabetes, arteriosclerosis or ischemia.
- People with hematological dyscrasias which affect thrombosis such as paroxysmal cold hemoglobinuria, cryoglobulinemia, sickle-cell disease or serum cold agglutinins.
- People who are in the acute stages of inflammatory phlebitis in the affected region.

(5) OPERATING THE SYSTEM


1. Push the lid release button to open the cooler lid (1), add water to fill-line indicated on label within the cooler (2). ENSURE WATER IS ABOVE THE FILL-LINE. Add as much ice as you want but do not exceed the ice max mark (3). After ice is added close the lid (4).


2. Apply the selected cuff (knee cuff as example). DO NOT APPLY DIRECTLY TO SKIN. Add an isolation layer between the cuff and the skin (pyjamas, jogging bottoms).

3. Connect cuff to cooler using connector hose. For best results, place the cooler at the same level as the cuff. The height difference between the two should not be more than 40cm.


4. Connect the AC adaptor to the cooler (12), then plug the AC adaptor (13).

5. Press the button for about 3-5 seconds, the screen should light up and show *Welcome*.

6. Press the or button to choose the mode.

Mode	Water circulation cycles	Time on (min)	Time off (min)
1	Low (ON: 4s / OFF: 12s)	20	20
2		40	40
3		60	60
4	Mid (ON: 8s / OFF: 7s)	20	20
5		40	40
6		60	60
7	High (ON: 12s / OFF: 7s)	20	20
8		40	40
9		60	60


Mode: There are 9 modes. Shows which program you have selected.

Press: There are 3 levels: low, medium, high. Represents the time of cycles of water circulation in the cuff during operation.

Time on: Time of active cycles of water circulation.

Time off: Pause time (no cycles of water circulation).

7. Press the button to start or pause the cooler. When cooling you can press the button to pause.
8. Press the button for about 3-5 seconds to turn off the cooler.
9. After use, to drain the cuff, disconnect the cuff from the cooler. Connect to the cuff the drain gadget. Once the drain gadget is connected squeeze the cuff to push the water out through the drain gadget (14). To drain the cooler, remove the lid and pour out the ice and water. Alternatively, you can use the drain gadget (15).



(6) CLEANING AND STORAGE

- Clean the surface of the cooler by wiping using a soft and damp cloth.
- Carefully dry the cuff with a dry cloth before storage to prevent from moisture.
- The product should be stored in a cool, dry, well ventilated and clean environment without corrosive substance.

(7) CUFFS

Shoulder Cuff

Ordering code: DA03020002


360 Knee Cuff

Ordering code: DA03020003



Ankle Cuff

Ordering code: DA03020004


Universal Cuff

Ordering code: DA03020005


Knee Cuff

Ordering code: DA03020010

Back Cuff

Ordering code: DA03020011

(8) ACCESSORIES

<p>Portable bag Ordering code: DA03020006</p>	
<p>Power bank Ordering code: DA03020007</p>	
<p>Battery case Ordering code: DA03020008</p>	
<p>Drain gadget Ordering code: DA03020009</p>	
<p>AC Adaptor Ordering codes: DA03020012 (UK) DA03020013 (Europe) DA03020014 (US) DA03020015 (AU/NZ)</p>	

(9) WARRANTY

The Kinetec® guarantee is strictly limited to the free replacement or factory repair of the part or parts that are recognized to be defective.

KINETEC SAS guarantees Dr Aktive® Cryo Pro for 2 years against any manufacturing defect, from the date of purchase by the consumer.

KINETEC SAS is solely accredited to judge the application of the guarantee for its devices.

The guarantee cannot be applied if the device has been used abnormally or was used in conditions of use other than those contained in our user manual.

The guarantee shall also not apply in the event of deterioration or accident resulting from negligence, compliance or maintenance fault resulting from a transformation of the equipment or an attempt to repair equipment.

(10) DISPOSAL AND RECYCLING

- **Packaging:** The packaging must be separated from plastic and paper/cardboard components and placed in specific recycling locations.
- **Cooler:** It contains electronic components, cables, plastic, steel and aluminium parts. When the device is no longer functional, dismantle and separate into groups of materials and place them in approved recycling units or return the device to KINETEC SAS for destruction. Or contact local authorities to establish the appropriate method of disposing of parts and accessories that are potentially dangerous to the environment.
- **Cuff and hose:** Put it in specific recycling locations.

(11) TROUBLESHOOTING

Cooler will not turn on	<ol style="list-style-type: none"> 1. Please make sure that AC adaptor is correctly inserted into the electrical outlet. 2. Please make sure that the power plug is correctly inserted into the cooler. 3. If you operate with battery, please make sure that the power bank is fully charged.
There is water on the surface of the compress pad	<ol style="list-style-type: none"> 1. Please make sure whether it is condensed water. 2. If the cuff leaks, stop using it immediately and contact the distributor for repair.
Loud noise when cold treatment starts	<ol style="list-style-type: none"> 1. Running for 2-3 minutes, the noise will be reduced. 2. If the noise is always very loud, please stop using it and contact the distributor for maintenance.

(12) SPECIFICATIONS

Weight of the cooler: 1.99kg (empty)

Capacity: maximum 3.6L

Acoustic pressure: < 50dBA

Power supply voltage: Input: 100-240V ~ 50/60Hz
Output: 12V = 2A

Electric type and class: BF Type / Class II

Storage/transport conditions:

Ambient temperature: -25°C to +70°C / -13 to 158°F

Relative humidity: .. from 20% to 90% without condensation

Atmospheric pressure: from 700hPa to 1060hPa

Operating conditions:

Ambient temperature: +5°C to +40°C / -13 to 158°F

Relative humidity: .. from 15% to 90% without condensation

Atmospheric pressure: from 700hPa to 1060hPa

(13) SYMBOLS USED

	Follow the instructions for use		Manufacturer
	CE marking of conformity		Medical device
	Keep away from sunlight and heat source		Device reference
	Keep dry		Device serial number
	Class II device (double insulation)		Manufacture batch number
	BF TYPE device (protection from electric shocks)		Manufacture date
	Alternating current		Direct current
	Unique Device Identifier information		Unique Device Identifier carrier
	Contains electric and electronic components; do not dispose of with household waste.		

(14) INFORMATION ABOUT ELECTROMAGNETIC COMPATIBILITY (EMC)

All the information given below comes from standard requirements to which electromedical devices are subject as defined by the IEC60601-1-2 standard.

The device radio-frequency emissions are very low and are not therefore likely to cause interference with electronic equipment installed nearby (radios, computers, phones, etc.). Nevertheless, users will make sure that possible electromagnetic interferences do not create an additional risk, such as from radio-frequency emitters or other electronic devices.

The device is designed to withstand foreseeable interferences from electrostatic discharges, magnetic fields from mains power supplies or radio-frequency emitters. Nevertheless, some types of mobile telecommunications devices such as mobile phones may interfere with the medical device. The separation distances recommended in this section must therefore imperatively be respected.

In this section, you will find the information required to guarantee the installation and commissioning of your medical device in the best conditions in terms of electromagnetic compatibility.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Dr Aktive® Cryo Pro, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

• Electromagnetic emissions

The device is designed for use in the electromagnetic environment described in the table below.

The user and installer must make sure the device is used in the environment described below.

Emissions tests	Compliance	Electromagnetic environment - Remarks
Electromagnetic radiation interference (radiated emissions) CISPR 11	Group 1	The device is suitable for use in a home healthcare environment and in a professional healthcare establishment environment.
Interference voltage at the supply terminals (conducted emissions) CISPR 11	Class B	
Harmonic current emissions IEC 61000-3-2	Class A	
Voltage variations, voltage fluctuations and flicker IEC 61000-3-3	Compliant	

• Magnetic and electromagnetic immunity

The device is designed for use in the magnetic and electromagnetic environment described in the table below.

Users and installers must guarantee the compliance of the electromagnetic environment.

Immunity tests	Test level as per IEC60601	Compliance level	Electromagnetic environment - Remarks
Electrostatic discharges (ESD) IEC 61000-4-2	±8 kV in contact ±2,4,8,15kV in the air	±8 kV in contact ±2,4,8,15kV in the air	The device is suitable for use in a home healthcare environment and in a professional healthcare establishment environment.
Fast electric transients / bursts IEC 61000-4-4	±2 kV for electricity supply lines	±2 kV for electricity supply lines	
Shock waves IEC 61000-4-5	±0.5, 1 kV in differential mode	±0.5, 1 kV in differential mode	The quality of the electricity supply grid must be equivalent to that of a home healthcare environment and a professional healthcare establishment environment.
	±0.5, 1, 2 kV in common mode	±0.5, 1, 2 kV in common mode	
Magnetic field at the assigned industrial frequency IEC 61000-4-8	30A/m	30A/m	The magnetic field intensity must be of the same level as that of a home healthcare environment and a professional healthcare establishment environment.
Voltage dips, short power cuts and voltage variations IEC 61000-4-11	0% U_T for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% U_T for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	The quality of the electricity supply grid must be equivalent to that of a home healthcare environment and a professional healthcare establishment environment. If the use of the device requires continued use during mains power outages, it is recommended to supply the medical system using a separate power source (UPS, etc.)
	0% U_T for 1 cycle	0% U_T for 1 cycle	
	70% U_T for 25 cycles at 50Hz for 30 cycles at 60Hz single phased at 0°	70% U_T for 25 cycles at 50Hz for 30 cycles at 60Hz single phased at 0°	
Voltage interruptions IEC 61000-4-11	0% U_T for 250 cycles at 50Hz for 300 cycles at 60Hz	0% U_T for 250 cycles at 50Hz for 300 cycles at 60Hz	

REMARK: U_T corresponds to the alternating network voltage before application of the test level.

• Electromagnetic immunity, portable RF equipment

The device is designed for use in the magnetic and electromagnetic environment described in the table below.

Users and installers must guarantee the compliance of the electromagnetic environment.

Immunity tests	Test level	Compliance level	Electromagnetic environment - Remarks
Radiated radio-frequency electromagnetic fields IEC 61000-4-3	3V/m from 80MHz to 2.7GHz 80%MA at 1kHz	3V/m from 80MHz to 2.7GHz 80%MA at 1kHz	The device is suitable for use in a home healthcare environment and in a professional healthcare establishment environment.
	10V/m from 80MHz to 2.7GHz 80% MA at 1kHz	10V/m from 80MHz to 2.7GHz 80% MA at 1kHz	
Proximity fields emitted by wireless radio-frequency communications devices IEC 61000-4-3	9 V/m 710MHz,745MHz, 780MHz,5240MHz, 5550MHz,5785MHz	9 V/m 710MHz,745MHz, 780MHz,5240MHz, 5550MHz,5785MHz	The device is suitable for use in a home healthcare environment and in a professional healthcare establishment environment.
	27 V/m 385MHz	27 V/m 385MHz	
	28 V/m 450MHz,810MHz, 870MHz,930MHz, 1720MHz,1845MHz, 1970MHz,2450MHz	28 V/m 450MHz,810MHz, 870MHz,930MHz, 1720MHz,1845MHz, 1970MHz,2450MHz	
Conducted interference induced by radio-frequency fields IEC 61000-4-6	3V from 150kHz to 80MHz	3V from 150kHz to 80MHz	The device is suitable for use in a home healthcare environment and in a professional healthcare establishment environment.
	6V in the ISM band and band included between 0.15MHz and 80MHz, including the amateur radio band 80%MA at 1kHz	6V in the ISM band and band included between 0.15MHz and 80MHz, including the amateur radio band 80%MA at 1kHz	

The electromagnetic field intensity of static radio-frequency emitters as determined by an electromagnetic environment measurement (a), must be less than the conformity level for each frequency range.

Interference can occur near equipment identified by the following symbol. 

REMARK: These specifications may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and persons.

(a) The electromagnetic field intensity of static radio-frequency emitters such as mobile phone base stations (cell / cordless), mobile radios, amateur radios, AF/FM radio emissions and TV emissions cannot be accurately determined by the theory.

To assess the electromagnetic environment due to the static radio-frequency emitters, an electromagnetic environment measurement must be carried out. If the measured intensity of the radio-frequency field in the immediate product use environment exceeds the radio-frequency conformity level specified above, the product performances must be tested to check that they are compliant with specifications. If abnormal performances are found, additional measurements may be needed, such as reorienting or moving the product.

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Kinetec SAS
Kinetec Medical Products UK Ltd (UK RP)
Kinetec USA, Inc.